

All ceramic restoration

COMPOSITION OF DENTAL PORCELAIN

1-The feldspars 70-90%

Mixtures of potassium aluminosilicate and sodium aluminosilicate.

The feldspar fuses when it melts, forming a glass matrix (give translucency)

2-The quartz (silica) 11-18%

remains unchanged during the firing process and acts as a strengthening agent

3-Kaolin 1- 10%


is a hydrated aluminum silicate ($\text{Al}_2\text{O}_3 \cdot 2\text{SiO}_2 \cdot 2\text{H}_2\text{O}$) and acts as a binder, increasing the ability to mould the unfired porcelain. , it also increase opacity

4-Metal oxides


which provide the wide variety of colors of the porcelain; for example, oxides of iron act as a brown pigment, copper as a green pigment, titanium as a yellowish-brown pigment, and cobalt imparts a blue color

Fabrication procedure

- The manufacturer mixes the components, adds additional metal oxides, fuses them and then quenches the molten mass in water. The resultant product is known as a *frit*, and the process is known as *fritting*.
- A consequence of the rapid cooling is that large internal stresses build up in the glass, resulting in extensive cracking. This material can be ground very easily to produce a fine powder for use by the dental technician



Sintering
→



Requirements of dental porcelain

- Low fusing temperature
- Retains its form during firing
- Resistance to devitrification
- High viscosity

Jacket crown construction

Die preparation



Platinum foil matrix


● **Function**

- 1- support of porcelain
- 2- even radiation of heat during firing

● **Properties of platinum matrix**

- 1-Inert during firing
- 2-Melting temperature higher than porcelain
- 3-Supplied in soft sheet which is easily adapted on prepared tooth

Platinum foil construction





Manipulation of dental porcelain

Application



Condensation


Purpose of condensation


- 1-remove excess water
- 2-adapt paste to the required form
- 3-decrease porosity
- 4-increase strength

Methods of condensation

- 1-Vibration
- 2- Capillary attraction
- 3-Spatulation

Porcelain build up





Reinforced ceramic core systems

- support is provided by another ceramic material

The addition of alumina to the feldspathic glass increases the strength of the core