



INVESTING



Presented by

Dr. Mohammed Hosny

Lecturer of Crowns and Bridges
Faculty of Dental Medicine
Al-Azher University

☺ The process of **casting** consists of surrounding the wax pattern with a mold made of heat resistant **investment** material, eliminating the wax by heating, and then introducing molten metal into the mold through a channel called **the sprue**.

DENTAL CASTING INVESTMENTS

They are composed of:

- **Refractory material**..... *Silica*, responsible for the thermal expansion of the mold during wax elimination
- **Binder**..... *gypsum, or magnesium oxide and an ammonium phosphate compound.*
- **Modifiers**..... accelerators, retarders and reducing agents *e.g. (carbon).*

Ideal properties of investment materials

- Able to control of expansion
- Produce smooth castings
- Accurate surface reproduction
- Chemical stability at high casting temperatures
- Adequate strength to resist casting forces
- Sufficient porosity to allow gas escape
- Easy recovery of the casting

Types of investment materials

- Gypsum-bonded
- Phosphate-bonded
- Silicate-bonded



Gypsum-bonded investments

- Gypsum binder, with cristobalite or quartz refractory material
- Not chemically stable at temperatures above 1000°
- More porous
- Mixed with distilled water
- Produce smooth surface casting
- Used for casting of conventional type II, III, and IV gold alloys

Phosphate-bonded investments

- Binder: magnesium oxide and ammonium phosphate compound.
- Silica refractory material
- Mixed with colloidal silica
- Higher strength
- Higher expansion due to the special liquid
- Stable at high temperatures; suitable for metal ceramic alloys

Phosphate-bonded investments

- Carbon-containing materials are used for high gold or palladium alloys
- Carbon-free investments are used for base metal alloys
- Lower porosity
- Rougher surface

Expansion

Setting expansion

Occurs during setting of the material as a result of crystal growth


Hygroscopic expansion


Extension of setting expansion by adding water to the setting investment by:-

- Wet liner.
- Submerging the ring in a water bath at 37°C for 1 hour immediately after investing.


Thermal expansion

Solid-state phase transformations of silica. More controllable



A

THERMAL EXPANSION
Cristobalite investment

**B**

How to increase expansion

- Use of two ring liners
- Prolonged spatulation
- Storage in 100% humidity
- Hygroscopic technique
- Lower water-powder ratio in gypsum-bonded and increase special liquid in phosphate-bonded


Investing techniques

■ **Single investing**

1. Brush technique
2. Vacuum technique

■ **Double investing**

Application of wetting agent (surfactant)






Mixing of investment






Brush technique




Vacuum technique








Double investing



Separation of crucible former



Wax elimination (burnout)

Aims:

- Complete elimination of wax to obtain an empty mold
- Obtain adequate thermal expansion
- Elevates the temperature of the investment to prepare it for casting

Burnout temperatures

- Gypsum-bonded: 650-700°C
- Phosphate-bonded: 850-900°C



Thank You!

