

METALS USED FOR CASTING

Dr. Ahmad El=Kouedi

REQUIREMENTS OF CASTING ALLOYS

- ❑ High strength and hardness
- ❑ Tarnish and corrosion resistance
- ❑ Acceptable fluidity and castability
- ❑ Can be soldered or welded
- ❑ Can add porcelain if needed
- ❑ Biocompatible with soft and hard tissue
- ❑ Easily finished and polished
- ❑ Can be burnished

Classification of alloys

- According to noble metal content
- -high noble $\geq 60\%$, AuPd, PtAuPdAg
- -noble 25% with no requirement of Au, AgPd
- -base metal less than 25 noble, NiCr, CoCr
- -titanium alloys

Classification of alloys

- According to mechanical properties(gold alloy)
- -type I: soft, used for small four walled cavities
- -type II: medium, used for complex restoration
- -type III: hard, used in Cr & Br
- -type IV: extra hard, used with RPD

Comparison between gold and base metal alloys

Gold alloy

- Yellow colour
- 18.3gm/cm³
- 935-1150 C
- 1.6% casting shrinkage
- Low modulus of elasticity
- Ductile and can be burnished

Base metal alloy

- Silver colour
- 7-9gm/cm³
- 1150-1400 C
- 2.4 %casting shrinkage
- High modulus of elasticity
- Inductile and can not be burnished

SPRUIING

- It's the process of adding a sprue former to the wax pattern before investing and casting.
- The purpose of the sprue former is to create a channel through which the molten alloy may enter the mold during casting.

Requirements of the sprue former

- Allows the escape of molten wax
- Allows molten alloy to flow into the mold
- Alloy in sprue must remain molten for a longer time than the alloy in the mold (reservoir)

Types of sprues

- According to material:
- 1-metallic sprue former:
 - - rust proof to prevent contamination
 - - hollow or solid
 - - easily removed but if the solid is used and it is forgotten then no casting failure will occur


Types of sprues

- 2- plastic sprue former:
- - requires higher melting temperature than wax, thus residual plastic may remain causing some casting failures.
- Are rigid so do not distort like wax

Types of sprues

- 3- wax sprue former:
- - most widely used, has same melting temp of wax pattern therefore the errors possible are reduced
- Comes in several ready made diameters

Sprue formers




Sprue former diameter

- - as a general rule the diameter of the sprue former should be the same thickness or larger than the thickest part of the wax pattern
- -2.5mm for molars/CM
- -2mm for premolars/partial coverage
- -too large will distort wax pattern
- -too small will not allow molten alloy to flow, unless pressure is applied



- Sprue former is attached to the thickest non functional cusp
- For anterior restorations at the mid incisal area
- Point of attachment should be flared to avoid constriction and turbulence
- The length of the sprue former will depend on the length of the ring and the type of the investment material used.
- Gypsum investment 6-8 mm
- Phosphate bonded 3-4mm
- Accessory or venting sprues can be used to allow casting of thin section

- 
- Direct sprueing uses with gold alloy
 - Indirect sprueing uses base metal alloy to prevent oxides from contaminating the cast
 - Indirect sprueing can use feeder, horizontal and oblique sprues





Crucible former

- It is usually made of rubber can be made from wood, metal, plastic
- It creates a funnel at the end of the sprue to allow entry of the molten alloy

Crucible former



Crucible former



Crucible former

